



HCPCS
A7005



COMPARE TO	
PARI	RESPIRONICS
022F81	HS870

Medel® JetNeb Plus

Item #	Product Description	UOM
18260	JetNeb Plus	12/cs

- The Medel® JetNeb Plus breath enhanced reusable nebulizer utilizes a dual valve system improving clinical efficacy and treatment efficiency.
- Through superior product engineering, the Medel® JetNeb Plus uses an inspiratory valve that delivers more medication while the patient breathes in, and an expiratory valve that wastes less when the patient breathes out.
- Optimized average particle size of 3.4 µm for targeted deep lung therapy.
- Up to 92.6% of medication is aerosolized into fine particle mist.

specifications	
Liter Flow	5.1 L / min
Output	0.45 ml
Output Range	0.02 ml per min
Nebulization Rate	0.3 ml per min
Average Particle Size	3.4 µm
Fraction of fine Particles	92.6%
Maximum Filling Volume	8 ml
Minimum Filling Volume	2 ml