



18301G
870 Regulator
For use with: M2 through E Cylinders



18303G
540 Regulator
For use with: M60, M90, MM, H, J, T Cylinders



18301GM
870 Mini Regulator
For use with: M2 through E Cylinders



18303GM
540 Mini Regulator
For use with: M60, M90, MM, H, J, T Cylinders

Drive Regulators

Item #	Product Description	Settings	Outlet	Connection	Liter Flow Increments	UOM
18300GN	4 lpm Pediatric Regulator	0-4 lpm	Barb	CGA 870	0.03, 0.06, 0.12, 0.25, 0.5, 0.75, 1, 1.5, 2, 3, 4	1/bx
18301G	8 lpm Regulator	0-8 lpm	Barb	CGA 870	0.12, 0.5, 1, 1.5, 2, 2.5, 3, 4, 6, 8	1/bx
18301GM	8 lpm Mini Regulator	0-8 lpm	Barb	CGA 870	0.12, 0.5, 1, 1.5, 2, 2.5, 3, 4, 6, 8	1/bx
18302GN	15 lpm Regulator	0-15 lpm	Barb	CGA 870	0.12, 0.5, 1, 1.5, 2, 3, 4, 6, 8, 10, 15	1/bx
18302GMN	15 lpm Mini Regulator	0-15 lpm	Barb	CGA 870	0.12, 0.5, 1, 1.5, 2, 3, 4, 6, 8, 10, 15	1/bx
18302GDEL	15 lpm Elongated Regulator for Humidifier	0-15 lpm	DISS	CGA 870	0.5, 1, 2, 3, 4, 5, 6, 8, 10, 12, 15	1/bx
18303G	8 lpm Regulator	0-8 lpm	DISS	CGA 540	0.12, 0.5, 1, 1.5, 2, 2.5, 3, 4, 6, 8	1/bx
18303GM	8 lpm Mini Regulator	0-8 lpm	DISS	CGA 540	0.12, 0.5, 1, 1.5, 2, 2.5, 3, 4, 6, 8	1/bx
18304GN	15 lpm Regulator	0-15 lpm	DISS	CGA 540	0.5, 1, 2, 3, 4, 5, 6, 8, 10, 12, 15	1/bx
18304GMN	15 lpm Mini Regulator	0-15 lpm	DISS	CGA 540	0.5, 1, 2, 3, 4, 5, 6, 8, 10, 12, 15	1/bx
18307G	4 lpm Pediatric Regulator	0-4 lpm	DISS	CGA 540	0.03, 0.06, 0.12, 0.25, 0.5, 0.75, 1, 1.5, 2, 3, 4	1/bx

- Drive regulators feature a protective contents gauge and a built-in pressure relief valve for safety.
- Mini regulators offer the same function as standard regulators while facilitating an easier fit into oxygen carry bags.
- Easy to operate single selector knob.

- Lightweight and durable design.
- Meet or exceed accuracy standards for ASTM, American National and CGA.
- Inlet Pressure: 300 - 3,000 PSI.
- Outlet pressure: 50 PSI.
- UL Approved.

specifications

Warranty	3 Year
----------	--------

accessories

	UOM
SEAL100 Replacement Seal-Washer for Oxygen Regulators; Made from Brass and Viton	10/bag