



**18302GDEL**  
**870 Regulator**  
 For use with: M2 through  
 E sized cylinders



**18303GM**  
**540 Regulator**  
 For use with: M60,  
 M90, H, J, T Cylinders

## Elongated Regulator For Humidifier

### 870 REGULATOR

- 1P** 18302GDEL 870 Style with DISS Outlet, 1/bx  
 Elongated to Attach Humidifier  
 Bottle Without Elbow Attachment.

### 540 REGULATORS

- 1P** 18303G Nut Connector w/Diss Outlet, Green, 1/bx  
 Liter Flow Increments: .5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 7, 8
- 1P** 18304GN Nut Connector w/Diss Outlet, Green,  
 (Has replaced 18304G), 1/bx  
 Liter Flow Increments: .5, 1, 2, 3, 4, 5, 6, 8, 10, 12, 15
- 1P** 18307G Nut Connector w/Diss Outlet, Pediatric, Green,  
 Liter Flow Increments: .03, .06, .12, .25, .5, .75, 1,  
 1.5, 2, 3, 4

### 540 MINI REGULATORS

- 1P** 18303GM Nut Connector w/Barb Outlet, Green, 1/bx  
 Liter Flow Increments: .5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 7, 8
- 1P** 18304GMN Nut Connector w/Barb Outlet, Green,  
 (Has replaced 18304GMN), 1/bx  
 Liter Flow Increments: .5, 1, 2, 3, 4, 5, 6, 8, 10, 12, 15

### FEATURES:

- New click style flow control.
- Lightweight uni-body design.
- Meets or exceeds accuracy standards for ASTM, American National and CGA.
- UL Approved.



## Replacement O-Rings for Oxygen Regulators

**SEAL100** Made From Brass and Viton Rubber  
 Packed, 10/bag