



18710



18711



NOW
COMES WITH
BLACK
RUBBER
PROTECTIVE
BOOT!



18705



18700

Health-OX Finger Tip Pulse Oximeter

18710

Retail, 1/bx

18710B

Protective Boot for 18710, 1/ea

FEATURES:

- Four way display.
- Ideal for individuals that exercise.
- Portable, non-invasive spot check oxygen saturation and heart rate.
- Two AAA batteries can be continuously operated for 30 hours.
- Operation is simple and efficient.
- Nail polish tolerant.
- 2 year limited warranty.

Medi- OX

18711

1/bx

18711B

Protective Boot for 18711, 1/ea

FEATURES:

- Four way display.
- Portable, non-invasive spot check oxygen saturation and heart rate.
- Two AAA batteries can be continuously operated for 30 hours.
- Operation is simple and efficient.
- Nail polish tolerant.
- 2 year limited warranty.

Clip Style Fingertip Pulse Oximeter

18705

1/bx

18705RB

Protective Boot for 18705, 1/ea

FEATURES:

- Two way display.
- Portable non-invasive, spot-check, oxygen saturation of arterial hemoglobin (SpO2) and pulse rate of adult and pediatric patient at home, and hospital (including clinical use in an internist/surgery, anesthesia, intensive care and etc). Not for continuously monitoring.
- Uses AAA batteries.

Fingertip Pulse Oximeter

18700

Rubber Boot Style, 1/bx

18700Battery

3.6V Lithium Battery

1/2 AA Size, 1/ea

FEATURES:

- Two way display.
- A spot check monitor for blood oxygen saturation and heart rate using reflectance pulse oximetry (RPO).
- Long-lasting battery - over 1,000 hours.
- Highly resistant to many forms of motion.
- Completely nail polish tolerant, no matter how dark the color.
- 2 Year Limited Warranty.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

respiratory

drive

To order call toll free: 877.224.0946

WWW.DRIVEMEDICAL.COM