



## Handy-Ox

Item #	Product Description	UOM
18720	Handy-Ox Digital Color Pulse Oximeter	1/ea

- The Handy-Ox Pulse Oximeter is a compact, convenient, high quality, portable oximeter with programmable pulse rate and SpO<sub>2</sub> alarm settings.
- Includes a lithium battery, AC adapter, an adult fingertip probe (pediatric probe available, sold separately), easy to use desktop software, and a USB cord to permit recorded data (up to 24 hours) upload or real time data transfer to a desktop computer.
- Two way, multi-view color thin film transistor (TFT) liquid crystal display.
- 2-digit SpO<sub>2</sub> and 3-digit heart rate numeric output with bar graph option.
- Included lithium battery has a life of no less than 500 charges.
- Audio and visual low battery alarm.

### specifications

Warranty	Limited 2 Year
System Requirements:	
Operating System	Microsoft Windows CP / Vista / 7
Processor	Pentium IV 1.8G or above
Memory	256MB or above
Motherboard	Intel chipset recommended
Disk Drive	40GB or above

### accessories

18720-ADPROBE	Reusable Adult Fingertip Pulse Oximeter Probe	1/ea
18720-PDPROBE	Reusable Pediatric Fingertip Pulse Oximeter Probe	1/ea