



**EASY HEIGHT
ADJUSTMENT
WITH SELF
THREADING
KNOB!**

10218FRD-1



199



**PARTS
DIAGRAMS
ON PAGE 66**



161



**EASY HEIGHT
ADJUSTMENT
WITH SELF
THREADING
KNOB!**

Winnie Lite Supreme/Go Lite 3 Wheel Aluminum Rollators, 7.5" Casters with Loop Locks

199 Tan Plaid, 1/cs,
10218FBL-1 Flame Blue, (Has replaced 196), 1/cs
10218FRD-1 Flame Red, (Has replaced 195), 1/cs

FEATURES:


- Lightweight solid 7.5" wheels for indoor or outdoor use.
- Comes standard with carry pouch.
- Special loop lock made of internal aluminum casting operates easily and ensures safety.
- Handles are adjustable in height.
- Adjustable brakes.
- Optional wire basket available (Item #806).
- Optional basket ring available (item #160).
- Easy, one hand folding.
- Limited Lifetime Warranty
Lifetime Warranty on brake cable.




Tan Plaid
Bag



Tan
Frame



Flame
Red
Frame



Flame
Blue
Frame

Winnie Deluxe 3 Wheel Steel Rollator, 7.5" Casters with Loop Locks

161 Chrome, (Has replaced 162), 1/cs

FEATURES:

- Lightweight solid 7.5" wheels for indoor or outdoor use.
- Comes standard with carry pouch.
- Adjustable handle height.
- Special loop lock made of internal aluminum casting operates easily and ensures safety.
- Optional wire basket available (Item #806).
- Optional basket ring available (item #160).
- Easy, one hand folding.
- Limited Lifetime Warranty.
Lifetime Warranty on brake cable.

SPECIFICATIONS

Item #	Handle Height (highest)	Handle Height (lowest)	Unit Width	Unit Length	Seat (depth)	Seat (width)	Seat (height)	Product Weight	Carton Shipping Weight	Weight Limit
WINNIE LITE SUPREME	38"	32"	26.5"	24"	NA	NA	NA	12 lbs.	16 lbs.	300 lbs.
WINNIE DELUXE	38"	32"	26.5"	24"	NA	NA	NA	12 lbs.	16 lbs.	300 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.