

TRIDENT

NEW!

Available
1 Pack

HCPCS
CODE:
K0823

THIS CHAIR IS THE "PERFECT" RENTAL CHAIR!

- ♦ Very durable design, high performance
- ♦ Easy to service
- ♦ Easy access to parts and components
- ♦ No shroud to be damaged
- ♦ Electronics can be accessed quickly
- ♦ Low maintenance cost for refurbishing

Trident, Front Wheel Drive Power Chair

2850-18 18" Seat, 1/cs
 2850-20 20" Seat, 1/cs

FEATURES:

- Durable, reliable performance, easy to service.
- Adjustable height, width and angle, padded arms.
- Seat is offered with 18 x 18 and 20 x 18 seat sizes.
- Semi reclining back, 20.5 height with headrest.
- Sits on contact points, seat plate can be adjusted forward or backward for patient positioning, very stable.
- Easy access to clutch to engage or disengage motor control.
- PG controller can be moved to either side for positioning.
- No shroud, easy open access to all components for quick servicing.
- Quick release pins for easy seat removal.
- Front wheel drive system, great power and navigation capability.
- Footplate height and angle adjustable.
- Weight Capacity 300 lbs.



SPECIFICATIONS

ITEM #	2850-18, 2850-20
HCPCS CODE GROUP	K0823
DIMENSION (L X W X H)	42" X 24" X 44.5"
GROSS WEIGHT (WITH/WITHOUT BATTERY)	186 LBS./133 LBS,
BATTERY WEIGHT	53 LBS.
SEAT WEIGHT	51 LBS.
SEAT DIMENSIONS (D X W)	18" X 18"/18" X 20"
BACK CASTER	8" SOLID
DRIVE WHEEL	14" SOLID
BATTERY (INCLUDED)	36 AMP X 2 U1
CHARGER	4A OFF-BOARD
MAX SPEED	4.0 MPH
CRUISING RANGE	12.4~18.6 MILES
CLIMBING ANGLE	6~8°
TURNING RADIUS	32.2"
GROUND CLEARANCE	2"
CONTROLLER	PG VR2 50A
WEIGHT CAPACITY	300 LBS.
WARRANTY	LIFETIME ON FRAME, ONE YEAR ON ELECTRONIC CONTROLLER AND DRIVE TRAIN COMPONENTS, SIX MONTHS ON BATTERY

powered mobility

To order call toll free: 877.224.0946

drive