

NEW!



**REPLACEMENT
FILTER**

ASK ABOUT OUR PRIVATE LABEL PROGRAM!

drive

SportNeb™2 Compressor Nebulizer

1P 3050-2 1/ea

HCPCS code: E0570

FEATURES:

- Powerful piston pump is efficient, dependable and quiet
- Compact design includes a convenient carry handle for this lightweight, 3.53 lbs. unit
- Built-in cord storage
- Includes 2 filters and disposable neb
- 5 Year Limited Warranty

ACCESSORIES:

18031
18254
Neb Kit 500
Mask 001a
Mask 001p
22311
3050-2-Filter

Universal Neb Carry Bag
Reusable Neb Kit
Disposable Neb Kit
Aerosol Mask - Adult
Aerosol Mask - Pediatric
Dragon Mask - Pediatric
Filter

SPECIFICATIONS	3050-2
METHOD OF NEBULIZATION:	Compressor
PARTICLE SIZE:	< 5 µm
UNIT DIMENSIONS:	4.1 " (H) x 5.7 " (W) x 7.9 " (D)
WEIGHT:	3.53 lbs.
MAXIMUM PRESSURE:	30 PSI
OPERATING PRESSURE:	10-14 PSI
LITER FLOW:	4-8 LPM
ELECTRICAL REQUIREMENTS:	115VAC/60HZ

1P = 1-Pack Available

respiratory

To order call toll free: 877.224.0946

drive