




The Industry's Most
Popular Innerspring
Mattress



Flex-Ease™ Firm Support Innerspring Mattress

Item #	Product Description	Product Dimensions	Product HCPCS	UOM
3637-1FE	Firm Innerspring Mattress	76" (L) x 36" (W) x 6" (H)	EO271	1/ea
3637-2FE	Firm Innerspring Mattress	80" (L) x 36" (W) x 6" (H)	EO271	1/ea
3637-3FE	Firm Innerspring Mattress	84" (L) x 36" (W) x 6" (H)	EO271	1/ea

- Economical innerspring mattress with two polyester fiber toppers.
- The unit has a 1" high-density, firm, combustion-modified fiber topper on each side.
- The Masongard vinyl waterproof cover is anti-microbial/anti-bacterial.
- The mattress is fire-retardant and conforms to CFR 16 part 1633.

specifications

Cover Material	Fluid-Proof Vinyl
Weight Capacity	250 lbs.
Warranty	1 Year