



Compare to:	<i>drive</i> <sup>™</sup>	Homecraft <sup>®</sup>
	Whisper	Deltis
Weight	20.5 lbs.	23.2
Lowest	2.3"	2.75
Reclines	50°	40°

**SPECIFICATIONS:**

Overall Length (with backrest in a reclined position)	32.6"
Overall Width (seat incl. side flaps)	27.5"
Width of Seat	13.7"
Depth of seat	16.5"
Overall Length of Seat	22.4"
Sitting Height Min.	2.3"
Sitting Height Max.	18"
Height of Backrest	25.1"
Width of Backrest	12.5" - 14.1"
Backrest Reclines	50°
Width of Baseplate	12.2"
Length of Baseplate	22"
Load Capacity Max.	300 lbs.
Total Weight (excl. hand control)	20.5 lbs.
Shipping Weight	25.3 lbs.

**Whisper Automatic Bath Lift**

- 1P** 477150252 Grey with White Cover, 1/bx
- 1P** 477150312 Grey with Blue Cover, 1/bx

**FEATURES:**

- The bath lift weighs only 20.5 lbs.
- Bath lift seat lowers to 2.3".
- Bath lift seat height reaches 18".
- The backrest reclines to 50 degrees at its lowest position.
- Completely tool free set up makes it easy to install and operate.
- The completely water tight hand control floats, making it easy to access the control buttons while bathing.
- The easy to use quick charge lithium-ion battery is accommodated in the hand control.
- Two piece construction makes the lift easy to assemble/disassemble for transport.
- When not in use, back easily folds down for space-saving storage.
- Padded seat and backrest provides added comfort.
- Hygienic cutaway and slim-line profile for easy handling.
- 4 quick release suction cups.
- Headrest available as an option.
- 3 Year Limited Warranty on the frame, hand control and charger.  
(Excluded from the warranty are wearable parts and suction cups.)

**1P = 1-Pack Available**