




The Bellavita is ultra modern,  
and the lightest bath lift on  
the market.



#### spare parts

460900403	Charger for Bellavita		1/bx
415252	Plug + Washer for Cover	White	1/ea
415312	Plug + Washer for Cover	Blue	1/ea
415432	Plug + Washer for Cover	Grey	1/ea
460900213	Suction Cup with clip	Grey 2.95" (75 mm)	1/ea
460900220	Set of Suction Cups with Clips	Grey 4 x 2.95" (75 mm)	4/pk
460900600	Hand Control Including Storage Battery		1/ea
460900225	Swiveling Sucker Set for Bellavita	Clear	4/pk



## Bellavita Auto Bath Lifter

#### Item #

477200252

#### UOM

1/cs

- The lightest bathlift on the market at only 20.5 lbs.
- The Bellavita lowers to 2.3". (Figure A)
- Seat height reaches 18.8" for deeper tubs.
- The backrest reclines to 50° at its lowest position.
- Completely tool-free set up makes it easy to install and operate. When not in use, back easily folds down for space-saving storage. (Figure B)

#### specifications

Overall Dimensions	32.6" (L) (with backrest in a reclined position) x 27.5" (W) (seat incl. side flaps)
Seat Dimensions	13.7" (W) x 16.5" (D) x 22.4" (L)
Sitting Height	2.3" - 18.8"
Backrest Dimensions	25.1" (H) x 12.5" - 14.1" (W)
Backrest Reclines	50°
Baseplate Dimensions	12.2" (W) x 22" (L)
Load Capacity Max.	300 lbs.
Total Weight (excl. hand control)	20.5 lbs.
Shipping Weight	25.3 lbs.
Warranty	Limited 3 Year

- Watertight hand control floats. (Figure C)
- The easy to use quick charge lithium-ion battery is accommodated in the hand control.
- Two-piece construction makes the lift easy to assemble/disassemble for transport.
- Padded seat and backrest provide added comfort.
- Washable hygienic covers are sold separately.
- Hygienic cutaway and slimline profile for easy handling.
- 4 quick-release suction cups for anti-slip bath base.
- Headrest available as an option.

#### accessories

412100200	Rotating & Transfer Aid	White	1/ea
410200252	Headrest - White	White	1/ea
410200312	Headrest - Blue	Blue	1/ea
410200432	Headrest - Gray	Gray	1/ea
460900205	Side Flap		1/ea