



**Horizontal & Vertical Cross-Cuts
Create 670 Independent
Pressure Redistributing Foam Cells**



Multi-Ply Shearcare500 Pressure Redistribution Foam Mattress

Item #	Product Description	Product Dimensions	Product HCPCS	UOM
500SC-1-FB	Single Layer/Multi-zoned Foam Mattress	76" (L) x 36" (W) x 6" (H)	E0184	1/ea
500SC-2-FB	Single Layer/Multi-zoned Foam Mattress	80" (L) x 36" (W) x 6" (H)	E0184	1/ea
500SC-3-FB	Single Layer/Multi-zoned Foam Mattress	84" (L) x 36" (W) x 6" (H)	E0184	1/ea
500SC-1-RR-FB	Single Layer/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out	76" (L) x 36" (W) x 6" (H)	E0184	1/ea
500SC-2-RR-FB	Single Layer/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out	80" (L) x 36" (W) x 6" (H)	E0184	1/ea
500SC-3-RR-FB	Single Layer/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out	84" (L) x 36" (W) x 6" (H)	E0184	1/ea

- A pressure redistribution foam mattress specifically created for facilities with limited budgets without sacrificing the clinical efficacy of more expensive surfaces.
- This three-zoned mattress is comprised of 670 independent cells.

specifications

Cover Material	Fluid-Resistant Nylon Top Waterproof Vinyl Bottom
Mattress Layers of Foam	1
Side Rails	Yes
Dedicated Heel Section	Yes
Fire Rating	CFR 16 Part 1633 Federal Fire Code
Weight Capacity	300 lbs.
Warranty	3 Year Limited Foam 1 Year Cover

- This design allows for maximum immersion and envelopment.
- Superior pressure redistribution is achieved by varying cell depths throughout the mattress zones.
- Unique SSPT Rail System™ provides solid-edge support, allowing easy patient transfer to and from bed.
- Designated head and foot sections prevent installation errors.
- Top Cover is a 2-way stretch nylon which reduces heat buildup and provides a low shear surface that is anti-microbial and water-resistant.
- Concealed zipper and a Barrier Stop Over-Flap prevents liquids from contaminating mattress core.
- Meets Federal Fire Code CFR 16 part 1633.