

LOOP LOCK
ENSURES
SAFETY!



Blue



Red



Green

Go-Lite Deluxe Aluminum Rollator, Padded Seat, 6" Casters with Loop Locks

10208PSBL-1	Blue, 1/cs
10208PSGR-1	Green, 1/cs
10208PSRD-1	Red, 1/cs

FEATURES:

- Special lever lock made of internal aluminum casting operates easily and ensures safety.
- 6" casters with soft-grip tires are ideal for indoor and outdoor use.
- Tool Free removable padded back rest with cover for comfort.
- Comes standard with carry pouch (item #804).
- Manufactured with wrap-around welds providing additional strength and support.
- Comfortable 12" x 12" padded seat.
- Handles are adjustable in height.
- Optional basket (item # 830) available.
- Limited Lifetime Warranty.

EASY HEIGHT
ADJUSTMENT
WITH SELF
THREADING
KNOB!

WEIGHS ONLY 12 LBS.

NOW COMES
STANDARD
WITH POUCH!

OPTIONAL BASKET
AVAILABLE
MODEL # 830

SEE PAGE 55



10208PSRD-1

NEW SEAT LOCK PREVENTS
ROLLATOR FROM
ACCIDENTALLY FOLDING.

SEE PAGE 52 & 53 FOR
REPLACEMENT SEATS.




Designer Winnie, Mimi Lite Deluxe Aluminum Rollator, Padded Seat, 6" Casters with Loop Locks

543	Flame Red Frame, 1/cs
544	Flame Blue Frame, 1/cs
548	Tan Frame, Tan Plaid Bag, 1/cs

FEATURES:

- 6" casters with soft-grip tires are ideal for indoor and outdoor use.
- Comes standard with pouch (item # 804).
- Tool Free removable padded back rest for comfort.
- Special loop lock made of internal aluminum casting operates easily and ensures safety.
- Handles are adjustable in height.
- Comfortable padded seat.
- Optional basket available (item # 830).
- Easy, one hand folding.
- Limited Lifetime Warranty.



**PARTS
DIAGRAMS
ON PAGE 61**

SEE PAGE 52 & 53 FOR
REPLACEMENT SEATS.

OPTIONAL BASKET
AVAILABLE
MODEL # 830

SEE PAGE 55

SPECIFICATIONS

Item #	Handle Height (highest)	Handle Height (lowest)	Unit Width	Unit Length	Seat (depth)	Seat (width)	Seat (height)	Product Weight	Carton Shipping Weight	Weight Limit
10208-1 SERIES	38"	32"	23"	24"	12"	12"	21.5"	12 lbs.	16 lbs.	300 lbs.
MIMI LITE DELUXE	38"	32"	23"	24"	12"	12"	21.5"	12 lbs.	16 lbs.	300 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

rollators

drive

To order call toll free: 877.224.0946

WWW.DRIVEMEDICAL.COM