



Gel-infused Foam Layer Keeps Skin Cool while Enhancing Patient Immersion & Envelopment for the Ultimate in Long Term Care Pressure Redistribution



Multi-Ply 6500 Dynamic Elite Pressure Redistribution Foam Mattress

Item #	Product Description	Product Dimensions	Product HCPCS	UOM
6500-DE-1-FB	Multi-layered/Multi-zoned Foam Mattress	76" (L) x 36" x (W) x 6" (H)	E0184	1/ea
6500-DE-2-FB	Multi-layered/Multi-zoned Foam Mattress	80" (L) x 36" x (W) x 6" (H)	E0184	1/ea
6500-DE-3-FB	Multi-layered/Multi-zoned Foam Mattress	84" (L) x 36" x (W) x 6" (H)	E0184	1/ea
6500-DE-1-RR-FB	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter Cut-out	76" (L) x 36" x (W) x 6" (H)	E0184	1/ea
6500-DE-2-RR-FB	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter Cut-out	80" (L) x 36" x (W) x 6" (H)	E0184	1/ea
6500-DE-3-RR-FB	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter Cut-out	84" (L) x 36" x (W) x 6" (H)	E0184	1/ea

- Gel-infused foam channels heat away from the body significantly faster than standard foam, keeping skin dry to minimize maceration.
- Superior pressure redistribution achieved by providing more support and less cradling at higher levels of compression, distributing pressure over a much larger surface area and eliminating static fatigue.
- Protects vulnerable heel area with the patented Lura-Quilt™ foam and a 25 degree heel slope redistributing pressure to the entire leg and away from the heel.
- 15 year non-prorated warranty.

specifications

Cover Material	Fluid-Resistant Nylon Top Waterproof Vinyl Bottom
Mattress Layers of Foam	4
Side Rails	Yes
Dedicated Heel Section	Yes
Fire Rating	CFR 16 Part 1633 Federal Fire Code
Weight Capacity	400 lbs.
Warranty	15 Year Limited Foam 1 Year Cover