



15-Year Warranty, Ultra-Soft, Aerated Foam LTC Mattress



Multi-Ply 6500 Global Pressure Redistribution Foam Mattress

Item #	Product Description	Product Dimensions	Product HCPCS	UOM
6500-GL-1-FB	Multi-Layered/Multi-zoned Foam Mattress	76" (L) x 36" (W) x 6" (H)	E0184	1/ea
6500-GL-2-FB	Multi-Layered/Multi-zoned Foam Mattress	80" (L) x 36" (W) x 6" (H)	E0184	1/ea
6500-GL-3-FB	Multi-Layered/Multi-zoned Foam Mattress	84" (L) x 36" (W) x 6" (H)	E0184	1/ea
6500-GL-1-RR-FB	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter Cut-out	76" (L) x 36" (W) x 6" (H)	E0184	1/ea
6500-GL-2-RR-FB	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter Cut-out	80" (L) x 36" (W) x 6" (H)	E0184	1/ea
6500-GL-3-RR-FB	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter Cut-out	84" (L) x 36" (W) x 6" (H)	E0184	1/ea

- Multi-layered pressure redistribution mattress that is soft enough for your lighter patients and durable enough for your heavier patients (450 lbs).
- The top layer of Lura-Quilt™ foam is aerated which allows air to continuously pass through reducing heat build up.
- The sloped heel section redistributes pressure from the heel to the less vulnerable calf.
- The heel section is precision die-cut to permit air to circulate under the patient's heel, which keeps the skin cool and dry.
- The Lura-Tex™ bottom layer is channel cut to allow continuous gatching of the head and foot section without compromising the integrity of the foam.
- Unique SSPT Rail System™ provides solid-edge support, and facilitates safe patient transfers.
- Designated head and foot sections prevent installation errors.
- Top Cover is a 2-way stretch nylon which reduces heat buildup and provides a low shear surface that is anti-microbial and water-resistant.
- Concealed zipper and a Barrier Stop Over-Flap prevents liquids from contaminating mattress core.
- Meets Federal Fire Code CFR 16 part 1633.

specifications

Cover Material	Fluid-Resistant Nylon Top Waterproof Vinyl Bottom
Mattress Layers of Foam	4
Side Rails	Yes
Dedicated Heel Section	Yes
Fire Rating	CFR 16 Part 1633 Federal Fire Code
Weight Capacity	450 lbs.
Warranty	15 Year Limited Foam 1 Year Cover