



3 Differential Pressure Zones on a Head-to-Toe Reversible Surface



Multi-Ply ShearCare700 Pressure Redistribution Foam Mattress

| Item # | Product Description | Product Dimensions | Product HCPCS | UOM |
|---------------|---|----------------------------|---------------|------|
| 700SC-1-FB | Multi-layered/Multi-zoned Foam Mattress | 76" (L) x 36" (W) x 6" (H) | E0184 | 1/ea |
| 700SC-2-FB | Multi-layered/Multi-zoned Foam Mattress | 80" (L) x 36" (W) x 6" (H) | E0184 | 1/ea |
| 700SC-3-FB | Multi-layered/Multi-zoned Foam Mattress | 84" (L) x 36" (W) x 6" (H) | E0184 | 1/ea |
| 700SC-1-RR-FB | Multi-layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out | 76" (L) x 36" (W) x 6" (H) | E0184 | 1/ea |
| 700SC-2-RR-FB | Multi-layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out | 80" (L) x 36" (W) x 6" (H) | E0184 | 1/ea |
| 700SC-3-RR-FB | Multi-layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out | 84" (L) x 36" (W) x 6" (H) | E0184 | 1/ea |

- Multi-layered pressure redistribution foam mattress that features three-zones with a specifically designed top layer of precision die-cut foam.
- The ShearCare™ 700 increases air circulation throughout the mattress, reducing heat buildup between the patient and sleeping surface.
- Mattress is fully reversible end-to-end, increasing product longevity.
- Bottom layer of high density foam provides a firm foundation.
- Unique SSPT Rail System™ provides solid-edge support, allowing easy patient transfer to and from bed.
- Top Cover is a 2-way stretch nylon which reduces heat buildup and provides a low shear surface that is anti-microbial and water-resistant
- Concealed zipper and a Barrier Stop Over-Flap prevents liquids from contaminating mattress core.
- Meets Federal Fire Code CFR 16 part 1633.

specifications

| | |
|-------------------------|---|
| Cover Material | Fluid-Resistant Nylon Top Waterproof Vinyl Bottom |
| Mattress Layers of Foam | 3 |
| Side Rails | Yes |
| Dedicated Heel Section | Yes |
| Fire Rating | CFR 16 Part 1633 Federal Fire Code |
| Weight Capacity | 300 lbs. |
| Warranty | 5 Year Limited Foam 1 Year Cover |