



**NOW SHIPPED IN
A SMALLER SHIP
CARTON!**

750NB

D-Lite, Aluminum Rollator, 8" Casters with Loop Locks and Tool Free Removable Wheels

750NB Blue, (Has replaced 750B
and 750NGC), 1/cs
750NR Red, (Has replaced 750R), 1/cs
750NBK Black, (Has replaced 750BK), 1/cs

FEATURES:

- Tool Free removable front and rear wheels.
- Tool Free removable contoured, padded back rest for comfort.
- Large 8" casters are ideal for indoor and outdoor use.
- Loop locks. Easy release without "snap-back".
- Large "mushroom" type button release on back.
- One piece release knob makes it easy to adjust handle height.
- Cut out folding padded seat.
- Comes with carry pouch.
- Ergonomic handles are easy to grip.
- Handles are adjustable in height.
- Easy, one hand folding.
- Limited Lifetime Warranty.
(Under License for Patent 7,066,484 and 7,073,801)




**ASSEMBLED WITH
TOOL FREE
REMOVABLE WHEELS**



**NEW SEAT LOCK PREVENTS
ROLLATOR FROM
ACCIDENTALLY FOLDING.**


**SEE PAGE 52 & 53 FOR
REPLACEMENT SEATS.**



Blue

Red

Black



**PARTS
DIAGRAMS
ON PAGE 64**



**EASY HEIGHT ADJUSTMENT
WITH SELF
THREADING KNOB.**



**WHEEL IS EASILY
REMOVED WITHOUT
TOOLS.**

SPECIFICATIONS

Item #	Handle Height Highest	Handle Height Lowest	Unit Width	Unit Length	Seat (depth)	Seat (width)	Seat (height)	Product Weight	Carton Shipping Weight	Weight Limit
750N Series	37.5"	31.5"	24"	27"	14"	13.5"	22"	19 lbs.	24 lbs.	300 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.