



Finally! An Affordable, Yet Clinically Effective Option in the Adjustable, Skin Protection Cushion Category



COMPARE TO

Compare to ROHO® High Profile® & Low Profile® Single Compartment Cushions

ROHO, High Profile and Low Profile are registered trademarks of ROHO, Inc.



Balanced Aire Adjustable Cushion

Item #	Product Description	Product Dimensions	Product HCPCS	UOM
8047-16	Adjustable Air Cushion	16" (W) x 16" (D) x 4" (H)	E2622	1/ea
8047-16-2	Adjustable Air Cushion	16" (W) x 16" (D) x 2" (H)	E2622	1/ea
8047-18	Adjustable Air Cushion	18" (W) x 16" (D) x 4" (H)	E2622	1/ea
8047-18-2	Adjustable Air Cushion	18" (W) x 16" (D) x 2" (H)	E2622	1/ea
8047-1818	Adjustable Air Cushion	18" (W) x 18" (D) x 4" (H)	E2622	1/ea
8047-1818-2	Adjustable Air Cushion	18" (W) x 18" (D) x 2" (H)	E2622	1/ea
8047-20	Adjustable Air Cushion	20" (W) x 16" (D) x 4" (H)	E2622	1/ea
8047-20-2	Adjustable Air Cushion	20" (W) x 16" (D) x 2" (H)	E2622	1/ea
8047-2018	Adjustable Air Cushion	20" (W) x 18" (D) x 4" (H)	E2622	1/ea
8047-2018-2	Adjustable Air Cushion	20" (W) x 18" (D) x 2" (H)	E2622	1/ea

- The Balanced Aire Adjustable Cushion is our premier pressure management cushion.
- Offering dynamic positioning, unmatched stability and an advanced level of pressure ulcer prevention.
- Flexible interconnected, independent air cells allow for deep immersion into the cushion.
- Optimal pressure redistribution helps increase blood flow necessary for healthy tissue maintenance.
- Water-resistant, low-shear stretch cover with adjustable locking straps and non-skid bottom to secure cushion to the wheelchair.
- Hand pump to adjust the cushions to individual patient comfort.
- Available in 4" and 2" adjustable air cell heights.

specifications

Cover Material	Fluid-Resistant Stretch Top Non-Skid Bottom
Weight Capacity	300 lbs.
Warranty	1 Year Limited