



## Convo-Gel™ Cushion

| Item #   | Product Description                                      | Product Dimensions         | Product HCPCS | UOM        |
|----------|----------------------------------------------------------|----------------------------|---------------|------------|
| 8049-U1  | 3" Gel/Convuluted Foam Cushion (Replaces Mason # 8049-1) | 16" (W) x 16" (D) x 3" (H) | E2603         | 1/ea, 4/cs |
| 8049-U2  | 3" Gel/Convuluted Foam Cushion (Replaces Mason # 8049-2) | 18" (W) x 16" (D) x 3" (H) | E2603         | 1/ea, 4/cs |
| 8049-U3  | 3" Gel/Convuluted Foam Cushion (Replaces Mason # 8049-3) | 20" (W) x 16" (D) x 3" (H) | E2603         | 1/ea, 4/cs |
| 8049-U4  | 3" Gel/Convuluted Foam Cushion                           | 20" (W) x 18" (D) x 3" (H) | E2603         | 1/ea, 4/cs |
| 8049-U5  | 3" Gel/Convuluted Foam Cushion                           | 22" (W) x 18" (D) x 3" (H) | E2604         | 1/ea, 4/cs |
| 8049-U6  | 3" Gel/Convuluted Foam Cushion                           | 22" (W) x 20" (D) x 3" (H) | E2604         | 1/ea, 4/cs |
| 8049-U7  | 3" Gel/Convuluted Foam Cushion                           | 24" (W) x 18" (D) x 3" (H) | E2604         | 1/ea, 4/cs |
| 8049-U8  | 3" Gel/Convuluted Foam Cushion                           | 24" (W) x 20" (D) x 3" (H) | E2604         | 1/ea, 4/cs |
| 8049-U9  | 3" Gel/Convuluted Foam Cushion                           | 26" (W) x 18" (D) x 3" (H) | E2604         | 1/ea, 4/cs |
| 8049-U10 | 3" Gel/Convuluted Foam Cushion                           | 26" (W) x 20" (D) x 3" (H) | E2604         | 1/ea, 4/cs |

- The Convo-Gel™ cushion is specifically designed to redistribute pressure while providing patient comfort.
- Cushion incorporates high-density convoluted foam over and under a dual-chamber gel bladder.
- Fluid-resistant, low-shear stretch nylon cover with tie straps to secure to wheelchair.

## specifications

|                 |                                          |
|-----------------|------------------------------------------|
| Cover Material  | Half P/U, Half Vinyl                     |
| Weight Capacity | Up to 20" 275 lbs.   22" and up 500 lbs. |
| Warranty        | 1 Year Limited                           |