



Pommel Wedge Cushion with Stretch Cover

Item #	Product Dimensions	UOM
8160-P	16" (W) x 16" (D) x 4"-2" (H)	1/ea, 4/cs
8161-P	18" (W) x 16" (D) x 4"-2" (H)	1/ea, 4/cs
8162-P	20" (W) x 16" (D) x 4"-2" (H)	1/ea, 4/cs
8163-P	22" (W) x 18" (D) x 4"-2" (H)	1/ea, 4/cs

- Wedge design helps prevent patient from sliding forward.
- Firm foundation reduces pressure caused by sling-seat hammockings.
- High-density foam increases stability and comfort.
- Foam pommel controls sliding by eliminating posterior pelvic tilt.

specifications

Cover Material	Fluid-Resistant Stretch Nylon, Half P/U, Half Vinyl
Weight Capacity	250 lbs.
Warranty	1 Year Limited

