

COUGAR

**QUICK RELEASE
REAR AXLE!**



**OPTIONAL
ADJUSTABLE
HEIGHT ARM
(Sold Separately)**



AK5ADJARMKIT



**BACK FOLDS FORWARD
TO REDUCE SIZE**

Cougar Wheelchair

Ultralight Aluminum Wheelchair (K0005)

FEATURES:

- Ultra light, weighs less than 27 lbs. (K0005).
- Forward folding back to reduce size.
- Tool free back height adjustable from 17"-20.5".
- Seat depth adjustable from 16"-18".
- Multi-position back can adjust 15°.
- Rigidity bar provides support.
- **Quick release rear wheels.**
- Semi-Pneumatic tires on a 24" spoke wheel provide a smooth, true roll.
- 3" of depth adjustment in the overall wheelbase.
- Rear axles are height adjustable from 15"-19".
- Front casters are angle adjustable and height adjustable in 5 positions.
- Arm rest can be flipped allowing 2" width adjustment.
- Nylon upholstery is breathable, lightweight and comfortable.
- Adjustable aluminum wheel locks.
- Swing-away footrests swing in and out.
- Lightweight, composite footplates with heel loops come standard.
- Limited Lifetime Warranty.

Due to product improvements and changes, all images in this catalog are subject to change.

**Important Notice: Item numbers that display an asterisk (*) reflect an item that will be discontinued or package change.*