



Deluxe Electronic Muscle Stimulator

AMS-4 1/bx

FEATURES:

- Intensity control is adjustable from 0-80 mA, 500 ohm load.
- Pulse rate is adjustable from 2-150 Hz.
- Dual isolated channels.
- Pulse width is selectable: 100, 200, or 300 micro-seconds.
- Contraction time is adjustable from 2-30 seconds.
- Relaxation time is adjustable from 1-30 seconds.
- Ramp time is adjustable from 1-6 seconds.
- Metal belt clip.
- Comes complete with carrying case, 2 lead wires, 9-volt battery, 1 package of 4 pre-gelled electrodes and operating manual.
- 3 Year Limited Warranty.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.