

WEIGHS ONLY 19 lbs.!



CONVENIENT CUP HOLDER CAN EASILY ATTACH TO CHAIR! SOLD SEPARATELY MODEL # STDS1040S

ATC19-BL



FRAMES



UPHOLSTERY



Aluminum Transport Chair

- ATC17-BK 17" Black Frame and Black Upholstery, (Has replaced AFW17-BK), 1/cs
- ATC17-BL 17" Blue Frame and Black Upholstery, (Has replaced AFW17-BL and FW17RW-BL), 1/cs
- ATC17-RD 17" Red Frame and Black Upholstery, (Has Replaced AFW17-RD and FW17RW-RD), 1/cs
- ATC17-SL 17" Silver Frame and Black Upholstery, (Has replaced AFW17-SL and FW17RW-SL), 1/cs
- ATC17-GR* 17" Green Frame and Black Upholstery, 1/cs
- ATC19-BK 19" Black Frame and Black Upholstery, (Has replaced AFW19-BK), 1/cs
- ATC19-BL 19" Blue Frame and Black Upholstery, (Has replaced AFW19-BL), 1/cs
- ATC19-RD 19" Red Frame and Black Upholstery, (Has replaced AFW19-RD), 1/cs
- ATC19-SL 19" Silver Frame and Black Upholstery, (Has replaced AFW19-SL and FW19RW-SL), 1/cs
- LELR-TF Swing-away Elevated Legrest, Tool Free, (Optional front riggings), 1pr/cs

FEATURES:

- Weighs only 19 lbs.
- Aluminum frame is lightweight and strong.
- Composite 8" wheels are lightweight and maintenance free.
- **Both 17" and 19" come with swing-away footrests.**
- Padded armrests provide additional comfort.
- Nylon upholstery is comfortable, lightweight and easy to clean.
- Aluminum casted back release hinge allows back to fold down for easy storage and transport.
- Standard rear wheel locks.
- Limited Lifetime Warranty.



SEE PAGE WC 278

SPECIFICATIONS

Item #	Overall Width (open)	Overall Width (closed)	Overall Height	Seat To-Floor Height	Back Height	Width Seat	Depth Seat	Overall Length w/ Rigging	Product Weight	Carton Shipping Weight	Weight Limit
ATC17	20"	9"	39.5"	19"	18"	17"	16"	33"	19 lbs.	24 lbs.	300 lbs.
ATC19	22"	9"	39.5"	19"	18"	19"	16"	33"	19 lbs.	24 lbs.	300 lbs.

Due to product improvements and changes, all images in this catalog are subject to change.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.