



**ASSEMBLES IN SECONDS!**



**ENVIRONMENT FRIENDLY PRODUCT**

Recycled materials were used to make this product.



## READY SET GO K.D. Deluxe, Aluminum, Bath Bench with Carry Bag

- 12105KDR-2 With Back, Knock Down, Retail, Comes with Carry Bag, 2/cs  
12105-4\* With Back, Knock Down, Individually Poly Bagged, Comes with Carry Bag, 4/cs  
12002T Replacement Tips for 12105, 4/bx

### FEATURES:

- Easy, tool free assembly of legs and back.
- Blow molded bench and back provides comfort and strength.
- Drainage holes in seat and back reduce slipping.
- Aluminum frame is lightweight, durable and corrosion proof.
- Angled legs with suction style tips provide additional stability.
- Retail model comes with a durable nylon carry bag for easy and convenient transport.
- Limited Lifetime Warranty.


**SUPPORT COLLAR PREVENTS LEG MOVEMENT!**

**RETAIL MODEL COMES WITH CONVENIENT CARRY BAG.**



**ASSEMBLED WITH TOOL FREE REMOVABLE BACK**

**LEGS AND BACK EASILY ASSEMBLED WITHOUT TOOLS!**



### Bath Bench Arms

12002A\* 1/pr

#### FEATURES:

- Easily attaches to removable back bath bench to provide additional stability. No tools required.
- For use with 12202-4 bath bench.



**SEE ACCESSORY PAGE 135 FOR UNIVERSAL BATHER POUCH ITEM # 10267**

### SPECIFICATIONS

| Item #     | Width (outside legs) | Width Seat | Depth (outside legs) | Depth Seat | Seat Height (highest) | Seat Height (lowest) | Product Weight | Carton Shipping Weight | Weight Limit |
|------------|----------------------|------------|----------------------|------------|-----------------------|----------------------|----------------|------------------------|--------------|
| 12105KDR-2 | 19"                  | 19.75"     | 18"                  | 12"        | 18"                   | 14"                  | 8.2 lbs.       | 19 lbs.                | 300 lbs.     |
| 12105-4    | 19"                  | 19.75"     | 18"                  | 12"        | 18"                   | 14"                  | 8.2 lbs.       | 35.8 lbs.              | 300 lbs.     |

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.

bath safety

To order call toll free: 877.224.0946

**drive**  
WWW.DRIVEMEDICAL.COM