



AquaJoy Premier Plus Reclining Bathlift

AquaJoy Saver Fixed-Back Bathlift

Item #	Product Description	UOM
BL100-DR	Reclining Back	1/ea

Item #	Product Description	UOM
BL200-DR	Fixed Back	1/ea

- Makes getting in and out of the bathtub safe and simple.
- Heavy duty seat and base provides maximum stability and support. (Figure A)
- Fits most style bathtubs.
- Fully waterproof, floating battery hand control with easy-to-use convex and concave operating buttons. (Figure B)
- Easy-release suction cups on base provide maximum stability in the tub.
- Comfortable, machine washable fabric covers.
- Safe. Will not lower unless there is sufficient power to raise the user back out of the tub.

- Makes getting in and out of the bathtub safe, simple and worry-free.
- Non-reclining, sturdy one-piece design maximizes legroom. (Figure A)
- Fits most style bathtubs.
- Fully waterproof, floating battery hand control with easy-to-use convex and concave operating buttons. (Figure B)
- Easy-release suction cups on base provides maximum stability in the tub.
- Safe. Will not lower unless there is sufficient power to raise the user back out of the tub.

specifications	
Seat Height Range	3.1" - 18.1"
Seat Width w/Flaps	27"
Recline	up to 40°
Powerlift and Recline	Yes
Weight	28.9 lbs.
Carton Shipping Weight	39.5 lbs.
Weight Capacity	375 lbs.
Warranty	Limited 3 Year

specifications	
Seat Height	2.7" - 18.3"
Seat Width w/Flap	28.5"
Weight	24.2 lbs.
Carton Shipping Weight	32 lbs.
Weight Capacity	350 lbs.
Warranty	Limited 3 Year