



Bobcat X4

Transportable Scooter

Item #	Seat	Color	UOM
BOBCATX4	16.5" Folding Seat	Red	1/ea

- Ultra-lightweight for superior portability.
- Quick and easy disassembly.
- Color-through shroud hides scuffs and scratches.
- Angle adjustable flip-up armrests for maximum comfort and convenient boarding and exiting.
- Easy finger-pull throttle control.
- Height-adjustable swivel seat with fold-down backrest.
- Flat-free, non-marking tires for worry-free travel.
- Large, metal carry basket.
- Easy-to-adjust tiller angle places the scooter's controls in the ideal position.
- Anti-tip wheels for additional safety.
- Easy freewheel operation.

specifications

Weight Capacity	265 lbs.	Charger	1.5A Offboard
Top Speed	4 mph	Brakes	Electromagnetic
Maximum Range	7.5 miles	Wheels	8" x 2" Flat-Free
Turning Radius	53"	Base Weight	59 lbs.
Climbing Angle	6°	Battery Weight (pair)	18 lbs.
Ground Clearance	2.5"	Seat Weight	16 lbs.
Dimensions (L x W)	42" x 19.25"	Total Weight	93 lbs.
Floor-to-Seat Height	17"-19"	Heaviest Piece Weight (Front Section)	38 lbs.
Seat Dimensions (W x D)	16.5" x 13.5"	Warranty on Frame	Lifetime
Controller	PG S-Drive 45A	Warranty on Electronic Controller/Drive Train Components	14 months
Motor	24V x 270W x 4700 rpm	Warranty on Batteries	6 months
Batteries (pair)	12V x 10AH		