



## Automatic Blood Pressure Monitor, Wrist Model

Item #	Product Description	UOM
BP2116	Automatic Blood Pressure Monitor, Wrist Model	1/ea

- Large easy-to-read digital display.
- 120 memory recall.
- Date & time display.
- Pulse Rate Indicator.
- Auto-off function.
- Lightweight & easy to use.
- Quiet inflation method.
- Wrist model.

### specifications

Measurement Method	Oscillometric
Measurement Range	Pressure: 0~300mmHg +/-3mmHg
Power Supply	2 x AAA Batteries
Auto Power Off	3 Minutes After Last Reading
Memory Capacity	120 Memories with Date & Time
Operation Environment	50F tp 104F, <85%RH
Warranty	2 Years

## Automatic Blood Pressure Monitor, Wrist Model

Item #	Product Description	UOM
BP2200	Automatic Blood Pressure Monitor, Wrist Model	1/ea

- WHO Classification.
- Irregular Heart Beat Detection.
- One Touch Operation.
- Large Easy-to-Read Screen.
- Date & Time Display.
- Pulse Rate Indicator.
- Travel Case.
- Auto-Off Function.
- Four User Memory Capable.

### specifications

Warranty	2 Years
----------	---------