

NEW!



CrawlAbout Crawl Trainer



CAB 2000 CrawlAbout, Small, Purple, 1/cs
CAB 3000 CrawlAbout, Medium, Blue, 1/cs

FEATURES:

- Height adjustable aluminum frame.
- Four easy glide pivoting casters.
- Height adjustable foam padded suspended harness. Adjust to user's size and ability to bear weight.
- Cushioned pad can support therapist or caregiver while they work with the user.
- Push handles can be attached as user's crawling skill progresses. The caregiver or therapist can steer from the rear to provide necessary assistance.
- Two year limited warranty.



ITEMS	OVERALL WIDTH	OVERALL HEIGHT	OVERALL DEPTH	HARNESS WIDTH	HARNESS DEPTH	WEIGHT	SHIPPING WEIGHT	WEIGHT LIMIT
CAB2000	23"	22"-25"	18"	12"	9"	10 lbs.	12 lbs.	50 lbs.
CAB3000	23"	26"-29"	23"	14"	11"	11 lbs.	13 lbs.	75 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

mobility
crawlabout

Wenzelite Re/hab

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM