



430

mobility
safety rollers
accessories

Vertical Handgrips

only available on posterior models

CE 1045 For use with all posterior Safety Rollers, 1pr/bx

Patients who cannot grasp a standard handlebar use vertical handgrips. The vertical handgrips extend, pronate, supinate and have a 45-degree angle. Each side is independently adjusted for the most comfortable hand placement. When both vertical handgrips are placed in the prone position, they can be used as a horizontal handlebar.



Pelvic Stabilizer

CE 1050 For use with all Wenzelite Safety Rollers, 1/bx

The pelvic stabilizer is adjustable in width and features a support belt. It is recommended for users who have an asymmetrical posture when standing or walking and, therefore, cannot stay centered in the walker.




Positioning Bar

CE 1052 For use with CE 1200 BK, 1/bx

CE 1053 For use with PE 1200, 1/bx

CE 1054 For use with PE TYKE 1200, 1/bx

The positioning bar is mounted on the handlebar of a Safety Roller and is easily moved for adjustable handlebar depth. It can be used with or without the pelvic stabilizer.



Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

WenzeliteRe/hab

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM