

Ankle Prompts

CE 1060 For use with all Wenzelite Safety Rollers, 1pr/bx

Maintains legs and/or feet in abduction, prevents crossing and controls stride length. Used for patients with spasticity, poor motor control, ataxia and excessive adduction.



Seat Harness

CE 1070S For use with all pediatric Safety Rollers & Nimbos, 1/bx

CE 1070L For use with all adult Safety Rollers & Nimbos, 1/bx

The Seat Harness moves with user's gait pattern allowing the user to shift weight from one leg to the other. The straps can be lengthened, shortened and positioned towards the front or back on the walker frame. It assists in weight bearing, properly positions user's pelvis and can be used as a safety sling. It features a seat pad that provides cushioned comfort.



Trunk Support

CE 1080 S for use with all pediatric Safety Rollers, 1/bx

CE 1080 L for use with all adult Safety Rollers, 1/bx

The trunk support wraps around the user and stabilizes the trunk. Used for patients with poor trunk control, ataxia, poor balance, spasticity or weakness.

FEATURES:

- Height and depth adjustable trunk support includes padded back and lateral pads. The lateral pads are width adjustable.
- Adjustable straps with front pad secures user for maximum support.
- Mounts directly on handlebar of both anterior and posterior Safety Rollers.



mobility
safety rollers
accessories

WenzeliteRe/hab

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.