

NEW! Nimbo

ACCESSORIES:

Pelvic Stabilizer

KA 1050 For use with All Nimbo Safety Rollers, 1/bx

The pelvic stabilizer is adjustable in width and features back & lateral pads and a support belt. It is recommended for users who have an asymmetrical posture when standing or walking and, therefore, cannot stay centered in the walker.



Nimbo Forearm Platform

KA 1035 FPS Forearm Platform, and Mounting Bracket, Small, (Has replaced KA1035, KA1035S), 1pr/bx

KA 1035 FPL Forearm Platform, and Mounting Bracket, Large, (Has replaced KA1035, KA1035S), 1pr/bx

- Constructed of lightweight aluminum.
- Height and depth adjustable.
- Handgrips extend forward or back, supinate, pronate and are angle adjustable.
- Individually adjustable
- Contoured armrest with a flexible edge.
- VELCRO® straps to prevent arms from slipping.
- Choice of mounting positions - inside or outside of the handlebar for lateral adjustment. Spacers are included for additional lateral range.



SPECIFICATIONS

Item #	Armrest Length	Height (from handlebar to armrest)
KA 1035 FPS	5"	3.5" - 4.5"
KA 1035 FPL	8"	5" - 7.5"

Seat Harness

CE 1070S Seat Harness, Small, 1/bx

CE 1070L Seat Harness, Large, 1/bx

The Seat Harness moves with user's gait pattern allowing the user to shift weight from one leg to the other. The straps can be lengthened, shortened and positioned towards the front or back on the walker frame. It assists in weight bearing, properly positions user's pelvis and can be used as a safety sling. It features a seat pad that provides cushioned comfort.



**RETROFITS
ON EXISTING
WALKERS!**



Seat

KA 1285 For use with foldable seat for KA 1200N, 1/ea

KA 2285 For use with foldable seat for KA 2200N, 1/ea

KA 3285 For use with foldable seat for KA 3200N, 1/ea

KA 4285 For use with foldable seat for KA 4200N, 1/ea

Seat mounts on standard Nimbo walker. They fold up for standing and walking and down for convenient seating.

