



Item #	Color	Seat	UOM
COBALTX23BL16FS	Blue	18" Folding Seat	1/cs
COBALTX23RD16FS	Red	18" Folding Seat	1/cs



HCPCS
K0823



Cobalt X23

Standard Power Wheelchair
Rear-Wheel Drive

- Full-size performance in a lightweight standard power wheelchair.
- Height-adjustable swivel seat features fold-down backrest.
- Armrests are padded, adjustable and removable.
- Adjustable-length controller mount.
- Optional swing-away controller arm (Item #AA4100) sold separately.
- Adjustable, flip-up footplate.
- Optional elevating leg rests (requires Item #AE2500 & Item #LK3JELR) sold separately.
- Positioning belt included.
- Armrest-mounted reflectors.
- Programmable controller.
- Flat-free, non-marking tires.
- Easy freewheel operation.
- Anti-tip wheels.

specifications

HCPCS	K0823	Batteries (pair)	12V x 21AH
Weight Capacity	250 lbs.	Charger	2A Offboard
Top Speed	4 mph	Brakes	Electromagnetic
Maximum Range	8 miles	Caster Wheels	8" x 2" Solid
Turning Radius	29"	Drive Wheels	9" x 3" Flat-free
Climbing Angle	6°	Base Weight	74 lbs.
Ground Clearance	2"	Battery Weight (pair)	28 lbs.
Dimensions	38.5" (L) x 24" (W)	Seat Weight	19 lbs.
Floor-to-Seat Height	20" - 23"	Total Weight	121 lbs.
Seat Dimensions	18" (W) x 15" (D)	Warranty on Frame	Lifetime
Controller	Dynamic Shark 40A	Warranty on Electronic Controller/Drive Train Components	14 months
Motors (2)	24V x 120W x 4700 rpm	Warranty on Batteries	6 months