

Star Posterior Gait Trainer

- SR 3000 Adult Star Posterior Gait Trainer, 1/cs
- SR 3100 Pediatric Star Posterior Gait Trainer, 1/cs
- SR 3200 Tyke Star Posterior Gait Trainer, 1/cs

Comet Anterior Gait Trainer

- CO 2000 Adult Comet Anterior Gait Trainer, 1/cs
- CO 2100 Pediatric Comet Anterior Gait Trainer, 1/cs
- CO 2200 Tyke Comet Anterior Gait Trainer, 1/cs

GAIT TRAINERS FOR PATIENTS WHO REQUIRE ANKLE PROMPTS TO PREVENT SCISSORING.

FEATURES:

- Ankle Prompts - prevents legs from scissoring.
- Seat Harness - supports patients who cannot stand completely on their own.
- Pelvic Stabilizer - supports and stabilizes the pelvis.
- Positioning Bar - provides adjustable depth.
- Vertical handgrips - provides alternate handgrip position, posterior models only
- Variable Resistance Tabs - controls the speed of the gait trainer and aligns an asymmetrical gait pattern.
- One-Directional Rear Wheels - allows only forward mobility, posterior models only.
- Swivel Wheel Locking Bracket - provides option to lock swivel wheels.
- Limited Lifetime Warranty on the frame.

ACCESSORIES:

- CE 1035 FP Adjustable Forearm Platforms and Mounting Brackets, 1pr/bx
- CE 1080 S Trunk support for pediatric & tyke anterior & posterior gait trainers, 1/bx
- CE 1080 L Trunk support for adult anterior & posterior gait trainers, 1/bx
- CE 1260 Extended uprights for the SR 3100, SR 3200, CO 2100 and CO 2200, 1pr/bx
- CE 1270 Extended uprights for the SR 3000 and CO 2000, 1pr/bx




Posterior



Anterior

ITEM #	ADJUSTABLE HEIGHT	INSIDE WIDTH BETWEEN HANDLEBAR	WIDTH	DEPTH OPENED AT BASE	WEIGHT LIMIT
SR 3000	29" - 36"	21"	26.5"	27"	400 lbs.
SR 3100	22" - 28"	15.5"	23.5"	25"	200 lbs.
SR 3200	17" - 20"	14"	22"	22"	100 lbs.
CO 2000	29" - 36"	21"	26.5"	27"	400 lbs.
CO 2100	22" - 28"	15.5"	23.5"	25"	200 lbs.
CO 2200	17" - 20"	14"	22"	22"	100 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.



mobility
gait trainers

Wenzel^{ite}Re/hab

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM