



HPCPS
A7037

CPAP-TUB 6

**Ask Us How to Private Label
Your CPAP Tubing**



HPCPS
A7037

CPAPTUB6TL

**Ask Us How to Private Label
Your CPAP Tubing**



Economy CPAP Tubing

Item #	Product Description	UOM
CPAP-TUB 6	Economy CPAP Tube, 6'	12/cs
RTLCAPTUB6	Retail Economy CPAP Tube, 6'	12/cs

- Durable and flexible 6' long CPAP tubing.
- Smooth bore universal tube.
- 22mm cuff connections are over-molded onto the tubing.
- Retail packaged available.

Trimline CPAP Tubing

Item #	Product Description	UOM
CPAPTUB6TL	Trimline CPAP Tube, 6'	12/cs
RTLCAPTUB6TL	Retail Trimline CPAP Tube, 6'	12/cs

- Premier 6', trim (15mm diameter) CPAP tubing.
- 15 mm diameter significantly reduces noise from PAP machine improving the patient therapy experience.
- 15 mm diameter also significantly reduces drag on mask, reducing leaks and improving patient comfort and compliance.
- 22mm cuff connections are over-molded onto the tubing.
- Retail packaged available.