

Back2Go Back Support



CAN BE MOUNTED TO ANY CHAIR!

A lightweight padded back support promotes improved posture and trunk control. It securely straps to chairs, wheelchairs and strollers.

CRB 2000 1/cs

FEATURES:

- Contoured back maximizes support and comfort.
- Two-strap system with buckle closure secures back support to an existing back surface. Straps can be stored in back pocket when not in use.
- Removable lateral supports are individually height and width adjustable.
- Zippered vinyl trunk harness in an H configuration for anterior chest support.
- 1 Year Limited Warranty



DIMENSIONS	CRB 2000
WIDTH	14"
HEIGHT	16"
DISTANCE BETWEEN LATERALS	6" - 9"
LATERALS	5" - 5"
SEAT TO LATERALS	4" - 7"



ACCESSORIES:

Headrest

CRS 8000 Headrest, 1/bx

Height adjustable headrest.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.