

CHROME SPORT



Chrome Sport Wheelchair Dual Axle (K0001/K0002/K0007)

FEATURES:

- **Chrome frame.**
- Padded armrests.
- Embossed vinyl upholstery is durable, lightweight, attractive and easy to clean.
- Urethane tires mounted on composite wheel provides durability and low maintenance.
- Comes with push to lock brakes.
- 8" front casters, adjustable in 3 positions.
- Carry pocket on backrest provides additional convenience.
- Comes with swing-away footrests or elevating legrests with tool free length adjustable front riggings.
- Plastic footplates.
- Color: black upholstery and chrome frame.
- 3 Year Limited Warranty.

**DURABLE
CHROME FINISH!**



Due to product improvements and changes, all images in this catalog are subject to change.

**Important Notice: Item numbers that display an asterisk (*) reflect an item that will be discontinued or package change.*