



D577-R

250 LB
WEIGHT CAPACITY

NEW & IMPROVED

Clinical Care Recliner

- D577-BR Clinical Care Recliner, Blue Ridge, 1/cs
- D577-J Clinical Care Recliner, Jade, 1/cs
- D577-R Clinical Care Recliner, Rosewood, 1/cs

AVAILABLE COLORS:



Rosewood Blue Ridge Jade



FEATURES:

- Ideal for use in dialysis, blood collection, respiratory and cardiac care.
- Chair can be positioned in 4 positions—upright, TV, Full Recline and Trendelenburg.
- Deep recline and Trendelenburg position can only be attained by an attendant.
- Gas cylinder controls deep recline and Trendelenburg position.
- Side panels "Pop-off" for easy cleaning and maintenance.
- Large, blow molded side trays includes recess for cup.
- 5" Casters (2 with locks) make it easy to move and secure recliner.
- Gas cylinder release handle for quick back release. Release handle can be operated from either side of the recliner.
- Wide, comfortable armrests are angled to provide support in any reclined position.
- Upholstered arm supports provide extra comfort in recline position.
- Neutral colored, sturdy steel frame is easy to clean and maintain.
- Comfortable built in headrest.
- Foam and vinyl meet California Flammability Standard and Boston Fire Department specifications.
- Limited Lifetime Warranty on Frame, 1 Year on Upholstery and Moving Parts.

SPECIFICATIONS

Item #	Overall Height	Overall Width	Overall Depth	Width Between Arms	Arm Height From Floor	Arm Height From Seat	Seat Height	Seat Width	Seat Depth	Carton Shipping Weight	Weight Limit
D577	47.5"	30.5"	38"	22"	26.5"	8"	22.5"	21.5"	19.5"	115 lbs.	250 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.