

Headrest

FC 8000 Adjustable Headrest for
FC 2000 & FC 4000, 1/bx

Height, depth and angle adjustable
headrest rotates and can be placed off center.

Headrest Wedges

WD 8000 Wedges, 1 pr/bx

Wedges can be placed anywhere on the
headrest for additional positioning.



Hip Guides

FC 8027 Hip Guides for FC 2000 & FC 4000, 1pr/bx

Hip Guides are width adjustable.

Mobility Base

FC 2029 Mobility Base for FC 2000, 1/bx

FC 4029 Mobility Base for FC 4000, 1/bx

Four wheeled mobility base safely transports child.

Leg Extensions

FC 8050 Leg extensions for FC 2000 & FC 4000, 1pr/bx

Set of four leg extensions for custom height.

Anti-tipper

FC 2030 Anti-tipper for FC 2000, 1 pr/bx

FC 4030 Anti-tipper for FC 4000, 1 pr/bx

Anti-tipper provides additional stability.



ITEM #	SEAT WIDTH	SEAT DEPTH	SEAT HEIGHT TO FLOOR	SEAT HEIGHT TO FOOTPLATE	BACK HEIGHT	LENGTH OF ARMREST	SEAT TO ARMREST HEIGHT	WEIGHT LIMIT
FC 2000	12"	9" - 14"	9" - 13"	6" - 12"	15"	10"	1.5" - 9.5"	100 lbs.
FC 4000	16"	13" - 17"	14" - 23"	6" - 22"	17"	10"	1.5" - 11.5"	200 lbs.
ITEM #	DISTANCE BETWEEN ARMRESTS	HIP GUIDES	DISTANCE BETWEEN HIP GUIDES	LATERAL	SEAT TO LATERALS	DISTANCE BETWEEN LATERALS	HEADREST	
FC 2000	9.5" - 14.5"	4" H X 9" D	8" - 13"	4" H X 5" D	1.5" - 6"	6" - 13"	ADDS 4.5" - 8" HEIGHT	
FC 4000	13.5" - 18.5"	4" H X 9" D	12" - 17"	4" H X 5" D	1.5" - 9"	10" - 17"	ADDS 4.5" - 8" HEIGHT	

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

positioning
seating accessories

Wenzel^{ite}Re/hab

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM