

HCPCS  
A7005



COMPARE TO			
DRIVE	HUDSON®	SALTER®	ROSCOE
18254	1901	8600	50016

COMPARE TO			
DRIVE	HUDSON®	SALTER®	ROSCOE
HUM 001	3260	7600	HUM-ROS6

## Reusable Nebulizer Kit

## Humidifier Bottle

Item #	Product Description	UOM
18254	Reusable Neb Kit with Mouthpiece, "T," and 7' Tubing	10/cs

- Anti-spill "T" design.
- Leak-proof and secure seal with simple 1/4 turn assembly design.
- Very low residual volume (0.7cc).
- High medication bowl capacity (10cc)
- Graduated volume (cc) markings on bowl.
- Nebulizes in upright positions or at angles up to 45°.
- Dishwasher safe.

Item #	For Use With	UOM
HUM 001	Oxygen Sources such as Concentrators, Liquid and Gas	50/cs

- Drive's multi-threaded lid ensures a secure fit every-time, preventing leaks and service calls.
- Audible back pressure alarm activates at 6 psi.
- No aluminum fitting in stem: it eliminates shredded aluminum diffuser so often found in lower quality models.
- Soft plastic diffuser reduces clogging due to mineral build-up.
- Translucent jar provides improved visibility of jar contents.
- Humidifier bottle lid is colored green/yellow, making it easy to tell patients on the phone how to change the bottle.
- Bottle has a 500 ml water capacity
- Bottle fill line is very easy to see.
- There are finger grips at the top of the bottle, making the lid easier to screw on/off for the patient.