


Nimbo

COMPARE & SAVE!

A less expensive alternative to the Kaye® posture walker!



Nimbo Lightweight Posterior Posture Walker

- KA 1200N Tyke, Goldenrod Yellow, 1/cs
- KA 2200N Junior, Lime Green, 1/cs
- KA 3200N Youth, Cornflower Blue, 1/cs
- KA 4200N Young Adult, Midnight Blue, 1/cs
- KA 5200N Adult, Flame Red, 1/cs **...NEW!**

NEWLY REDESIGNED!

FEATURES:

- Aluminum Frame in bright colors.
- Height adjustable.
- Soft rubber wheels adhere to any surface and will not allow the Nimbo to slide backwards.
- One directional rear wheels with aluminum ratchet and pin mechanism.
- Anti-Reverse override bracket disengages the one directional rear wheels and allows forward & reverse mobility.
- Front swivel wheels enhance maneuverability and ease of turning.
- 5" rear wheels (image not shown). **...NEW!**
- Locking mechanism easily switches front wheels from swivel to non-swivel. Each wheel is independently adjusted.
- The KA 2200N, KA 3200N, KA 4200N and KA 5200N can be used in both anterior and posterior position..
- Latex-free handgrips with flanged ends.
- Folds easily for convenient transporting.
- 3 year Limited Warranty.



Fixed Wheel

Swivel Wheel

*** For proper sizing, measure wrist to floor height.**

Many physical therapists prefer posterior (reverse) walkers for children with Cerebral Palsy and other neurological, orthopedic and developmental disorders.

ITEM #	HANDLEBAR HEIGHT RANGE	OVERALL WIDTH	WIDTH INSIDE HANDLE GRIP	DEPTH OPENED AT BASE	PRODUCT WEIGHT	CARTON SHIPPING WEIGHT	WEIGHT LIMIT
KA 1200N	14.5" - 18"	24"	13.5"	20"	9 lbs.	11 lbs.	75 lbs.
KA 2200N	19" - 25"	24"	13.5"	27"	10 lbs.	13 lbs.	100 lbs.
KA 3200N	23" - 30.5"	26"	15"	31"	13 lbs.	14 lbs.	150 lbs.
KA 4200N	28" - 36"	26.5"	16.5"	34"	14 lbs.	18 lbs.	200 lbs.
KA 5200N	36" - 41"	29"	19"	37"	16 lbs.	20 lbs.	275 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.


mobility
posterior posture walker

nimbo

Wenzel *ite Re/hab*

To order call toll free: 877.224.0946

drive



Nimbo

ACCESSORIES:

Pelvic Stabilizer

KA 1050 For use with all Nimbo Posture Walkers, 1/bx

The pelvic stabilizer is adjustable in width and features back & lateral pads and a support belt. It is recommended for users who have an asymmetrical posture when standing or walking and, therefore, cannot stay centered in the walker.



Forearm Platform

KA 1035 FPS Forearm Platform and Mounting bracket, Small, 1pr/bx

KA 1035 FPL Forearm Platform and Mounting bracket, Large, 1pr/bx

- Constructed of lightweight aluminum.
- Height and depth adjustable.
- Handgrips extend forward or back, supinate, pronate and are angle adjustable.
- Individually adjustable.
- Contoured armrest with a flexible edge.
- Velcro® straps to prevent arms from slipping.
- Choice of mounting positions - inside or outside of the handlebar for lateral adjustment. Spacers are included for additional lateral range.

ITEM #	ARMREST LENGTH	HEIGHT (FROM HANDLEBAR TO ARMREST)
KA 1035 FPS	5"	3.5" - 4.5"
KA 1035 FPL	8"	5" - 7.5"



Dial-A-Speed Tabs

KA 1400 For Use with All Nimbos, 1pr/bx

Dial-A-Speed tabs provide variable resistance on one or both rear wheels. They allow user to control walker speed and aligns asymmetrical gait patterns.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

Seat Harness

CE 1070S Seat Harness, Small, 1/bx

CE 1070L Seat Harness, Large, 1/bx

The Seat Harness provides cushioned comfort. It moves with user's gait pattern allowing the user to shift weight from one leg to the other. The straps can be lengthened, shortened and positioned towards the front or back on the walker frame. It assists in weight bearing, properly positions user's pelvis and can be used as a safety sling.



Seat

KA 1285 For use with KA 1200N, 1/ea

KA 2285 For use with KA 2200N, 1/ea

KA 3285 For use with KA 3200N, 1/ea

KA 4285 For use with KA 4200N, 1/ea

KA 5285 For use with KA 5200N, 1/ea

Solid seats fold up for standing and walking and flip down for convenient seating.



ITEM #	SEAT-TO-FLOOR HEIGHT RANGE
KA 1285	9" - 13"
KA 2285	11.5" - 17"
KA 3285	12" - 19.5"
KA 4285	14.5" - 22"
KA 5285	21.5" - 27.5"

**RETROFITS
ON EXISTING
WALKERS!**

mobility

posterior posture walker accessories

nimbo

Wenzelite Re/hab

To order call toll free: 877.224.0946

drive