

## Soft Seat Harness



Nimbo

HCPCS  
E1399

Item #	Product Description	UOM
CE 1070S	Small	1/bx
CE 1070L	Large	1/bx

- The Seat Harness provides cushioned comfort. It moves with user's gait pattern allowing the user to shift weight from one leg to another.
- Straps can be lengthened, shortened and positioned towards the front or back of the gait trainer's frame
- It assists in weight bearing, properly positions user's pelvis and can be used as a safety sling.

## Pouch



Nimbo

HCPCS  
E1399

Item #	Product Size	UOM
10258-1	Small	1/bx
10258L-1	Large	1/bx



- Nylon carry pouch features two pockets.

## Seat



Nimbo

HCPCS  
E1399

Item #	For Use With	Product Seat-to-Floor Height Range	UOM
KA 1285	KA 1200N	9" - 13"	1/bx
KA 2285	KA 2200N	11.5" - 17"	1/bx
KA 3285	KA 3200N	12" - 19.5"	1/bx
KA 4285	KA 4200N	14.5" - 22"	1/bx
KA 5285	KA 5200N	21.5" - 27.5"	1/bx

- Solid seats fold up for standing and walking and flip down for convenient seating.

## 9" Front Wheels



Nimbo

HCPCS  
E1399

Item #	For Use With	UOM
KA 8100	KA 2200N-KA 5200N	1 pr/bx

- Non-swivel front wheels with inner legs.
- Height adjustable.
- Easily used outdoors and over rough grounds.
- Adds 1" to overall height of walker.