



HCPCS
E1236



Pediatric Kanga TS Inclusive System

Item #	Product Size	Color	Product Seat-to Floor Ht	Product Overall Width	Product Seat Depth	Product Overall Depth w/ Riggings	Product Weight	Product Weight Capacity	UOM
KG 1000	Pediatric, 10"	Red	21" (Front) - 20" (Back)	18" (Open) - 12" (Closed)	11.5"-13.5"	37"	47 lbs.	200 lbs.	1/bx
KG 1200	Pediatric, 12"	Red	21" (Front) - 20" (Back)	20" (Open) - 12" (Closed)	11.5"-13.5"	37"	47 lbs.	200 lbs.	1/bx
KG 1400	Pediatric, 14"	Red	21" (Front) - 20" (Back)	22" (Open) - 12" (Closed)	11.5"-13.5"	37"	47 lbs.	200 lbs.	1/bx

- Durable lightweight aluminum frame.
- Built-in growth frame is width adjustable and can be configured for 10", 12" and 14". Seating system and positioning accessories must be changed accordingly.
- Height- and angle-adjustable flip-back armrests.
- Height-adjustable push handles.
- Tilt-in-Space up to 45 degrees.
- 6" of adjustable seat depth.

- Seat-to-back angle adjustable.
- Angled seat to prevent forward sliding.
- Swing-away, height-adjustable footrests with foot and ankle positioner.
- Angle-adjustable footplates.
- 20" rear semi-pneumatic quick-release wheels.
- Frame accommodates most custom rehab seating and positioning systems.
- Curb-assist levers.
- Height-adjustable anti-tippers.
- Crash tested and meets bus transit requirements.
- Folds for easy transport. (Figure A)

specifications

Push Handle Height	38.75" - 43.5"
Armrest Height	10" - 12"
Back Cane Height	18"
Back Angle Adjustment	5° - 25°
Carton Shipping Weight	53 lbs.
Warranty	Limited Lifetime

accessories

KG 8024	Kanga TS Upper Extremity Tray	1/bx
KG 8026	Kanga TS Canopy	1/bx