

Kanga TS

Inclusive System

Features a seating system that includes seat, five-point harness, depth adjustable abductor, hip guides, lateral supports and curved headrest.



Model #:

Sold 1 per box

Pediatric

	KG 1000	10", Red
	KG 1200	12", Red
	KG 1400	14", Red

Adult

	KG 1600	16", Black
	KG 1800	18", Black
	KG 2000	20", Black

200/300 lbs

Limited Lifetime

Weight Capacity:

Warranty:

Accessories see pg 28-29

Features

- Durable lightweight aluminum frame.
- **Built-in growth frame is width adjustable. Pediatric model can be configured for 10", 12" and 14", the adult model for 16", 18" and 20".** Seating system and positioning accessories must be changed accordingly.
- Height and angle adjustable flip-back armrests.
- Height adjustable push handles.
- Tilt-in-Space up to 45 degrees.
- 6" of adjustable seat depth.
- Seat to back angle adjustable.
- Angled seat to prevent forward sliding.
- Swing away, height adjustable footrests with foot and ankle positioner.
- Angle adjustable footplates.
- 20" rear semi-pneumatic quick release wheels.
- Frame accommodates most custom rehab seating and positioning systems.
- Curb assist levers.
- Height adjustable anti-tippers.
- Crash tested and meets bus transit requirements.
- Folds for easy transport.
- Limited Lifetime Warranty.

Specifications

Item #	Overall Width (open)	Overall Width (closed)	Front Seat-to Floor Ht	Back Seat-to Floor Ht	Seat Depth	Overall Depth w/Riggings	Product Weight	Carton Shipping Weight	Weight Limit
KG 1000	18"	12"	21"	20"	10"-16"	37"	47 lbs.	53 lbs.	200 lbs.
KG 1200	20"	12"	21"	20"	10"-16"	37"	47 lbs.	53 lbs.	200 lbs.
KG 1400	22"	12"	21"	20"	10"-16"	37"	47 lbs.	53 lbs.	200 lbs.
KG 1600	24"	12"	22"	21"	15"-21"	42"	52 lbs.	59 lbs.	300 lbs.
KG 1800	26"	12"	22"	21"	15"-21"	42"	52 lbs.	59 lbs.	300 lbs.
KG 2000	28"	12"	22"	21"	15"-21"	42"	52 lbs.	59 lbs.	300 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

 = 1-Pack Available