



KG P800



KG A800





HCPCS
E1234
(KG P800)

HCPCS
E1161
(KG A800)

Pediatric & Adult Kanga TS Folding Tilt-in-Space Wheelchair Frame

Offers Built-In Growth for Long Term Fit!

Item #	Product Size	Growth Frame	Color	Product Overall Width Open	Product Seat Depth	Product Weight w/ Wheels	Product Weight w/o Wheels	Product Weight Capacity	UOM
KG P800	Pediatric	Can be configured for 10", 12" & 14"	Red 	18" - 22"	13" - 15"	34 lbs.	26 lbs.	200 lbs.	1/bx
KG A800	Adult	Can be configured for 16", 18" & 20"	Black 	24" - 28"	17" - 19"	35 lbs.	27 lbs.	300 lbs.	1/bx

specifications

Overall Width (closed)	12"
Seat-to-Floor Height	19"
Back Cane Height	18"
Push Handle Height	38.75" - 43.5"
Armrest Height	10" - 12"
Overall Length with Footrest	37"
Overall Length Without Footrest	28.5"
Back Angle Adjustment	5° - 25°
Carton Shipping Weight (KG P800)	34.5 lbs.
Carton Shipping Weight (KG A800)	38 lbs.
Warranty	Limited Lifetime

- Durable lightweight aluminum frame.
- Height- and angle-adjustable flip-up armrests.
- Height-adjustable push handles.
- Tilt-in-Space up to 45°.
- Adjustable-width cross brace. (Figure A)
- 6" of seat depth.
- Seat-to-back angle adjustable.
- Curb-assist levers.
- 20" rear semi-pneumatic quick release wheels.
- Accommodates most custom rehab seating and positioning systems.
- Crash tested with Kanga TS seating and meets bus transit requirements.
- Folds for easy transport. (Figure B)