

VIPER



AVAILABLE IN 12" AND
14" SEAT WIDTHS!

NOW
COMES WITH
RIGGINGS IN
CARTON!

Viper
 Deluxe High Strength,
 Lightweight, Dual Axle (K0004)

FEATURES:

- VIPER weighs 32 lbs. (Excluding front rigging).
Easier to propel and transport.
- New frame style eliminates seat guides and allows for custom backs and accessories.
- Back height adjustable from 18" - 20".
- Carbon steel black powder coated frame double coated for an attractive, chip-proof, maintainable finish.
- Comes standard with "flip up" anti-tippers.
- Composite steel lined front receiver is angled so the arm will lock into place every time.
- The front mounted release lever can be flipped up or down depending on the patients needs.
- Durable, heavy gauge flame retardant nylon upholstery resists mildew and bacteria.
- Nylon, padded armrests provide added patient comfort.
- Flip-back and removable arms allow for easy transfer.
- Composite Mag-style wheels with aluminum hand rims are lightweight and maintenance free.
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability.
- Standard dual axle provides easy transition of seat height to hemi-level.
- Composite footplates with heel loops are durable and lightweight.
- Front caster forks are adjustable in 3 positions.
- Urethane rear tires and casters offer superior performance, smoother ride and are lightweight.
- Comes with push to lock wheel locks.
- Color: Black.
- Limited Lifetime Warranty.

*** NOT FEATURED ON 12" & 14" MODELS**

Due to product improvements and changes, all images in this catalog are subject to change.

**Important Notice: Item numbers that display an asterisk (*) reflect an item that will be discontinued or package change.*



RELEASE LEVER DOWN



RELEASE LEVER UP

FRONT MOUNTED RELEASE LEVER CAN
BE FLIPPED UP OR DOWN DEPENDING
ON THE PATIENT'S NEEDS.*

NEW &
IMPROVED
"SLOPE"
STYLE ARM

