

Luminator Posterior Gait Trainer

- LT 3000 Adult Luminator Posterior Gait Trainer, 1/cs
 LT 3100 Pediatric Luminator Posterior Gait Trainer, 1/cs
 LT 3200 Tyke Luminator Posterior Gait Trainer, 1/cs

Luminator Anterior Gait Trainer

- LT 2000 Adult Luminator Anterior Gait Trainer, 1/cs
 LT 2100 Pediatric Luminator Anterior Gait Trainer, 1/cs
 LT 2200 Tyke Luminator Anterior Gait Trainer, 1/cs

FEATURES:

- Seat Harness – supports patients who cannot stand completely on their own.
- Pelvic Stabilizer – supports and stabilizes the pelvis.
- Positioning Bar – provides adjustable depth.
- Vertical Handgrips – provides alternate handgrip position, posterior models only.
- Variable Resistance Tabs – controls the speed of the gait trainer and aligns an asymmetrical gait pattern.
- One-Directional Rear Wheels – allows only forward mobility, posterior models only.
- Swivel Wheel Locking Bracket – provides option to lock swivel wheels.
- Limited Lifetime Warranty on the frame.



Posterior



Anterior

ACCESSORIES:

- CE 1035 FP Adjustable forearm platforms and mounting brackets, 1pr/bx
 CE 1060 Set of ankle prompts, 1pr/bx
 CE 1080 S Trunk support for pediatric & tyke anterior & posterior gate trainers, 1/bx
 CE 1080 L Trunk support for adult anterior & posterior gait trainers, 1/bx
 CE 1260 Set of extended uprights for the LT 3100, LT 3200, LT 2100 and LT 2200, 1pr/bx
 CE 1270 Set of extended uprights for the LT 3000 and LT 2000, 1pr/bx

ITEM #	ADJUSTABLE HEIGHT	INSIDE WIDTH BETWEEN HANDLEBAR	WIDTH	DEPTH OPENED AT BASE	WEIGHT LIMIT
LT 3000	29" - 36"	21"	26.5"	27"	400 lbs.
LT 3100	22" - 28"	15.5"	23.5"	25"	200 lbs.
LT 3200	17" - 20"	14"	22"	22"	100 lbs.
LT 2000	29" - 36"	21"	26.5"	27"	400 lbs.
LT 2100	22" - 28"	15.5"	23.5"	25"	200 lbs.
LT 2200	17" - 20"	14"	22"	22"	100 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.