



Item #	Color	Seat	UOM
MEDALIST450BL22CS	Blue	22" Captain's Seat	1/cs
MEDALIST450RD22CS	Red	22" Captain's Seat	1/cs
MEDALIST450BL24CS	Blue	24" Captain's Seat	1/cs
MEDALIST450RD24CS	Red	24" Captain's Seat	1/cs



HCPCS

K0825

Medalist 450

Heavy-Duty Power Wheelchair
Rear-Wheel Drive

- Height-adjustable swivel seat features fold-down backrest.
- Four-post seating.
- Seat back reclines and features adjustable headrest.
- Adjustable-length controller mount.
- Armrests are padded, adjustable and removable.
- Available with 22" or 24" width Captain's Seat.
- Optional swing-away controller arm (Item #AA4100) sold separately.
- Adjustable, flip-up footplate.
- Positioning belt included.
- Optional elevating leg rests (requires Item# AE3000 & Item #HDELR) sold separately.
- Programmable controller.
- Flat-free, non-marking tires.
- Easy freewheel operation.
- Rear battery drawer allows easy access for service or replacement.
- Anti-tip wheels.

specifications

HCPCS	K0825	Batteries (pair)	12V x 55AH
Weight Capacity	450 lbs.	Charger	5A Offboard
Top Speed	4 mph	Brakes	Electromagnetic
Maximum Range	15 miles	Caster Wheels	9" x 2.5" Solid
Turning Radius	31"	Drive Wheels	14" x 3" Flat-free
Climbing Angle	6°	Base Weight	167 lbs.
Ground Clearance	3"	Seat Weight	55 lbs.
Dimensions	40" (L) x 24" (W)	Battery Weight (pair)	66 lbs.
Floor-to-Seat Height	20" - 23"	Total Weight	288 lbs.
Seat Dimensions	22" (W) x 20" (D) 24" (W) x 20" (D)	Warranty of Frame	Lifetime
Controller	Dynamic Shark 60A	Warranty on Electronic Controller/Drive Train Components	14 Months
Motors (2)	24V x 500W x 4600rpm	Warranty of Batteries	6 Months