





MICHAEL GRAVES
active living
COLLECTION™

MICHAEL GRAVES TUB RAILS



MG12036



MG12041



Clamp On Tub Rail

MG12036 White, Non-Adjustable, Parallel Design, Retail, 1/bx

MG12041 White, Non-Adjustable, Perpendicular Design, Retail, 1/bx

FEATURES:

- **MG Orange Activates.** Turn knob for easy adjustment.
- Steel construction for safety and stability.
- Designed for beauty as well as functionality.
- Easy tool-free assembly.
- Durable powder-coated finish.
- Fits standard size tub. Width adjusts from 3" - 7".
- No sharp edges makes bathing safer and more comfortable and won't scratch tub.
- Limited Lifetime Warranty.



MG12050SC



U.S. Patent # D566,250



Clamp On, Height Adjustable Tub Rail with Soft Cover, Soap and Shampoo Dish

MG12050SC White, Adjustable, Retail, 1/bx

FEATURES:

- **MG Orange Activates.** Turn knob for easy adjustment.
- Soft cover for added safety.
- Non-slip comfort grip makes bathing safer and more comfortable for everyone.
- Oversized handle provides additional stability, making bathing safe and secure.
- Easy tool-free assembly.
- Height adjusts from 14" - 17". Width adjusts from 3" - 7".
- Integrated dishes for soap and shampoo.
- No sharp edges makes bathing safer and more comfortable and won't scratch tub.
- Limited Lifetime Warranty.



SPECIFICATIONS

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

Item #	Adjustable Clamp (smallest)	Adjustable Clamp (widest)	Height (lowest to highest)	Product Weight	Carton Shipping	Weight Limit
MG12036/41	3"	7"	14" - 17"	6.4 lbs.	8.4 lbs.	300 lbs.
MG12050SC	3"	7"	14" - 17"	7.4 lbs.	9.4 lbs.	300 lbs.

bath safety

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM