



HCPCS
A7005



COMPARE TO			
DRIVE	HUDSON®	SALTER®	ROSCOE
18254	1901	8600	50016



Reusable Nebulizer Kit

Pirate Island Building Block Kit

Item #	Product Description	UOM
18254	Reusable Neb Kit with Mouthpiece, "T," and 7' Tubing	10/cs

- Anti-spill “T” design.
- Leak-proof and secure seal with simple 1/4 turn assembly design.
- Very low residual volume (0.7cc).
- High medication bowl capacity (10cc)
- Graduated volume (cc) markings on bowl.
- Nebulizes in upright positions or at angles up to 45°.
- Dishwasher safe.

Item #	Product Description	UOM
MQ0073	Pirate Island Building Block Kit for MQ7000 Compressor System - Compressor Not Included	1/ea