



Vital Sign Monitor

Item #	Product Description	UOM
MQ3600	Vital Sign Monitor	1/ea
MQ3600P	Vital Sign Monitor with Printer	1/ea
MQ3600T	Vital Sign Monitor with Thermometer	1/ea

- Compact, lightweight and portable.
- Easy Operation with low maintenance.
- SpO₂ Measurement with Memory.
- Pulse Rate Measurements.
- Automated Monitoring of Noninvasive Blood Pressure (NIBP).
- Data Storage and Adjustable Alarm Settings.

specifications

Clinic Mode	1000 Patients / 16000 Recordings
Single User Mode	1 - 16000
Power	Rechargeable Lithium ion Battery
Operating Time	8 hrs.
Charge Time	6 hrs.
Power Supply	100-24-v, 50/60Hz, OVA
Weight	5.51 lbs.
Warranty	Limited 1 Year

accessories

MQ3601	SpO ₂ ExtCable, A Series, 9Pin to 9Pin	1/ea
MQ3602	Digital SpO ₂ Ad ReuseSensor, 9 Pin	1/ea
MQ3603	Digital SpO ₂ Ped ReuseSensor, 9Pin	1/ea
MQ3609	Adult Cuff with Connector	1/ea
MQ3610	Large Adult Cuff with Connector	1/ea
MQ3612	Infant Cuff with Connector	1/ea